



## Workstation Specifications for Picis LYNX C/Point® and I/Point™ 6.2

### C/Point Workstations

| Minimum   | Recommended   |
|---|---|
| <ul style="list-style-type: none"> <li>• Intel® Pentium® III 750 MHz+</li> <li>• Microsoft Windows 2000 Professional SP4</li> <li>• Memory: 256 MB RAM</li> <li>• Available hard drive space: 500 MB</li> <li>• .NET Framework 2.0 sp1 or 3.5</li> <li>• Microsoft Internet Explorer 6.0</li> <li>• 15" Monitor</li> <li>• Mouse or other pointing device</li> <li>• 128 Kbps synchronous Internet speed or greater.</li> </ul> | <ul style="list-style-type: none"> <li>• Intel® Core™ 2 Quad 2.66 GHz+</li> <li>• Microsoft® Windows® XP SP2, SP3, Windows Vista®, Windows 7, Windows Server® 2003 or 2008</li> <li>• Memory: 4 GB RAM*</li> <li>• Available hard drive space: 1 GB</li> <li>• .NET Framework 2.0 sp1 or 3.5</li> <li>• Microsoft Internet Explorer 7.0 or later</li> <li>• 19" Monitor</li> <li>• Mouse or other pointing device</li> <li>• 640 Kbps synchronous Internet speed or greater.</li> </ul> |

\* Maximum RAM may be constrained by Operating System

### Printers

- HP LaserJet 4250 or higher (model to be determined based on facility annual census)
- 64 MB RAM (minimum)
- 100 Mbps NIC (unless connected locally)

## Delivering Results in Revenue Management™

Document Control #992025710A

Copyright 2010, LYNX Medical Systems, Inc. All rights reserved.

LYNX Medical Systems, logo, I/Point, and C/Point are trademarks or registered trademarks of LYNX Medical Systems, Inc. Intel and Pentium are trademarks of Intel Corporation in the U.S. and other countries. Microsoft, Windows, Windows Server, Windows Vista, Silverlight, and Internet Explorer are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries. Adobe and Reader are either registered trademarks or trademarks of Adobe Systems Incorporated in the United States and/or other countries.